



CUSTOM NOISE GENERATOR ORDER FORM


PHOENIX
Hearing Instruments
ABN 49 087 517 262

| | |
|------------|---------|
| Mr | _____ |
| Ms | _____ |
| Mrs | _____ |
| Miss | _____ |
| FIRST NAME | SURNAME |

Date: _____

Clinician: _____ Telephone: _____

| | |
|-------------|----------|
| Account To: | Ship To: |
| ABN: | |



PHOENIX
Hearing Instruments
Date Received _____

| |
|------------|
| Serial # L |
| Serial # R |

| |
|---------|
| Box |
| Invoice |

| HISTORY | Left | Right |
|---------------|--|--|
| Previous User | <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Years worn | <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Years worn |
| Model | _____ | _____ |
| Output | _____ | _____ |
| Gain | _____ | _____ |
| Vent | _____ | _____ |
| Years Used | _____ | _____ |

Tinnitus Retraining Therapy

L R **PHING** - Phoenix Noise Generator (Skeleton ITE style with Janssen vent)

L R **PHING** - Phoenix Noise Generator (Canal ITC style with trench vent)

“TINNITUS RETRAINING THERAPY (TRT) - is a treatment regime based directly on the Jastreboff model. It uses a precise and individual combination of sound therapy and teaching / demystification / learning about tinnitus and hyperacusis mechanisms, with directions about how to return to normal life without provoking symptoms.” (The Tinnitus and Hyperacusis Centre, London UK - Jonathan Hazell and Jacqueline Sheldrake)

The Phoenix noise generator (PHING) was developed to be used in support of a TRT based approach to Tinnitus and Hyperacusis management.

Because of the very nature of the condition resolution is often long term and sound therapy may be required for a protracted period in conjunction with the other aspects of therapy before benefit is observed. For this reason, devices are not accepted for credit return. Within the first 3 months, the usual provisions for remake for issues with comfort or poor fit will apply. A one year warranty applies and covers any defects, which under proper use, appear in the goods that are due to faulty materials, workmanship or design.

The PHING is designed to be as open fit as possible. ITC versions of the device are not recommended, however, some clients do find they are helpful and are able to overcome some comfort issues when sleeping with the smaller style of device.

Please ensure that a copy of the audiogramme is attached to your order. This way we can ensure there is adequate sensation level provided by the device.

N.B. This device is **NOT** a tinnitus masker and is not designed for use without appropriate adjunctive therapy.